English

Mouse Anti-Human Actinin Alpha 4 Monoclonal Antibody (Anti-ACTN4 mAb)

Catalog ID: DH0003

For in Vitro Diagnostic Use (IVD): USA, EU and Taiwan

For Research Use Only (RUO): Other countries

For Professional Users

Intended Use

Mouse Anti-Human Actinin, Alpha 4, clone 13G9 is intended for the semi-quantitative detection of Actinin, Alpha 4 protein in paraffin sections. The clinical interpretation of any positive staining or its absence should be complemented by morphological studies and histology studies with proper controls. Evaluations should be made within the context of the patient's clinical history and other diagnostic tests by a qualified pathologist.

Principles of Procedure

Immunohistochemical(IHC) staining techniques allow for visualization of tissue constituents localization under microscope through a two step process including antigen-primary antibody interactions and the detection of bound antibody by a chromogen.

Reagent Provided

1. 100 ug of monoclonal mouse antibody, clone 13G9, to Actinin Alpha 4, diluted in 1X PBS, pH 7.2.

Reagent Required but Not Provided

1. Blocking solution or Antibody Diluent (DAKO, Cat # S3022)
2. Immunodetection Kit- EnVision Detection Kit, Peroxidase/DAB,Rabbit/Mouse (DAKO, Cat # K5007)

Users are advised to use the reagents recommended by Abnova, otherwise the expected results may not be achievable.

Isotype

IgG2b

Immunogen

peptide, 'NQSYQYGPSSAGNGA'.

Specificity

Human Actinin, Alpha 4

Antibody Concentration

Refer to vial label for batch specific Ig concentration.

Recommendation on Working Condition

Suggest primary antibody concentration: 1.5ug/ml overnight at 4°C, and do not re-use. Heat mediated antigen retrieval using 0.01M citrate retrieval solution (pH 6.0) is recommended. These are guidelines only and the users should determine their own optimal condition.

Storage and Stability

Please refer to vial label for expiration date and store at -20°C. Do aliquot to avoid repeated freezing and thawing. Do not mix different lot of antibody into one and not use after expiration date. Storage condition other than those specified in the package insert, they must be verified by the user.

Specimen Preparation

Tissues specimens preserved for IHC staining by formalin fixation.

Staining interpretation

The cellular staining pattern for anti-Human Actinin, Alpha 4 is cytoplasm or nucleus.

Staining Procedure

The following protocol is used in ABNOVA. Customer can refer to the protocol from the detection system selected.

Deparaffinize sections and rehydrate using 1X PBS

Pre-treat the sample with the following procedure:

Place sample in 1X citrate buffer (pH 6.0) in pressure cooker under 125°C for 4min and under 90°C for 45min, cool sample subsequently. Wash sample with 1XPBS.

Step-by-step procedure:

1. Incubate sections in 3% H₂O₂ in 1X PBS at room temperature for 10 minutes and then wash the sections with 1X PBS.
2. Incubate sections in blocking solution for 10 minutes.
3. Add primary antibodies (1.5 ug/ml) and incubate the sections overnight at 4°C, wash sample with 1X PBS afterwards.
4. Incubate sections with labeled polymer⁽¹⁾ (HRP-conjugated 2 Ab) for 30 min followed by washing the sections

with 1X PBS.

5. Application of chromogen and substrate solution⁽²⁾ (DAB or other suitable peroxidase substrate). Wash sample thoroughly under running tap water.
6. Counter stain the samples in Mayer's hematoxylin.
7. Dehydrate and mount samples.

Note:

⁽¹⁾ ⁽²⁾ These are the components of the Immunodetection Kit- EnVision Detection Kit, Peroxidase/DAB, Rabbit/Mouse (DAKO, Cat # K5007)

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mdi Europa GmbH